




# ESRF Status Report

## Daniele de Sanctis

**On behalf of the rest of the team:**

Marcus Oscarsson, Antonia Beteva, Olof Svensson, Axel Bocciarelli, Jean Baptiste Florial (*EMBL*) and the SB-Group



- MXCuBE Web 3 on ID23-2 and MASSIF3
  - MX beamlines are running almost completely remotely
- MXCuBE Web 4 with mxcubecore deployed in production on
  - ID29 For SSX experiments
  - ID30A1 For Automatic data collection, Harvester and Plate integration
  - ID23-1 Mainly proprietary research now
  - ID30B currently commissioning
  - ID23-2 during winter shutdown
- BSXCuBE3 on BM29
- Since last MXCuBE meeting
  - Initial Pump&probe experiments on ID29
  - Implementation of in-situ (plates) experiments on MASSIF1 (N.Foos)
  - Deployment of GPhL workflows
  - Developments on ID23-1 to increase sample throughput

- MXCuBE Web 4
  - Based on **mxcube**core
  - Major Backend and frontend library **updates**
  - New **sample view**
  - **In situ** data collection
  - **Configurable display**
  - **Braggy** and heatmap integration
  - **SSX** modes for ID29
  - Adding **Methods** and procedures in three ways:
    - **Equipment view** - For not so often used or temporary instrumentation commands
    - **Beamline action** - For procedures that are frequently used and involves more than a simple command
    - **Queue entry / task** - For collecting data

MXCuBE-Web (osc)

Samples Data collection Equipment System log

Help Remote Sign out (idtest0)

Beamline Actions ▾

Energy: 12.4000 keV Resolution: 0.740 Å Transmission: 10.0 % Cryo: 200.00 K

Wavelength: 1.00 Å Detector: 10.1 mm Flux: 2.30e+12 ph/s

Camera 1 Camera 2 Detector Sample Changer Fast Shutter Safety shutter Ring Current

■ ■ ■ READY READY CLOSED CLOSED 200.0 mA

The display of available instrumentation is configurable in ui.yaml

To the left motor control and on the top “beamline setup”

MXCuBE-Web (os)

Beamline Actions ▾

Phase Control:  
Centring

Beam size:  
10

Omega:  
311.10 90.0 °

Kappa:  
11.0 0.1 °

Kappa Phi:  
22.0 0.1 °

Sample alignment:

⬆ ⬅ ⚙ ➡ ⬇

Show motors ▾

```
sample_view:  
  id: sample_view  
  components:  
    -  
      label: Omega  
      attribute: diffractometer.phi  
      role: omega  
      step: 90  
      precision: 2  
      suffix: °  
      # Marcus Oskarsson, 20 months ago * I  
    -  
      label: Kappa  
      attribute: diffractometer.kappa  
      role: kappa  
      step: 0.1  
      precision: 1  
      suffix: "°"
```

```
beamline_setup:  
  id: beamline_setup  
  components:  
    -  
      label: Beamstop  
      attribute: beamstop  
    -  
      label: Capillary  
      attribute: capillary  
    -  
      label: Fast Shutter  
      attribute: fast_shutter  
    -  
      label: Safety shutter  
      attribute: safety_shutter  
    -  
      label: Detector  
      attribute: detector  
    -  
      label: Energy  
      attribute: energy  
      step: 0.001  
      precision: 4  
      suffix: keV
```

## Equipment view - For less often or temporary instrumentation commands

MXCuBE-Web (osc)    Samples    Data collection    **Equipment**    System log    ? Help    Remote    Sign out (edtest)

Mockup **READY**

Content

Refresh    Scan all containers

+ Mockup

Power

PowerOn    PowerOff    Regulation On

Lid

Open Lid    Close Lid

Actions

Home    Dry    Soak

Recovery

Clear Memory    Reset Message    Back    Safe

Abort

Abort

**my\_fancy\_function**

my\_fancy\_function    Last result

Speed\*

Num Images\*

Exp Time\*

PhaseEnum\*  
An enumeration.

Run my\_fancy\_function

detector **READY**

diffractometer **READY**


abort    my\_fancy\_function

**Methods** are automatically added if they are “exported” with the export tag and the method is type hinted (at least with a return type)

```
<exports>["abort", "status", "my_fancy_function", "my_other_funny_function"]</exports>
```


**Beamline action - For procedures that are frequently used and involves more than a simple command**

**Configured as the `beamline_actions` of the Beamline hardware object**



```
<object class="BeamlineActionsMockup">
  <commands>[
    {"type": "controller", "name": "Centre Beam", "command": "HardwareObjects.m
    {"type": "controller", "name": "Quick Realign", "command": "HardwareObjects
    {"type": "controller", "name": "Anneal", "command": "HardwareObjects.mockup
    {"type": "annotated", "command": "HardwareObjects.mockup.BeamlineActionsMoc
    {"type": "annotated", "command": "HardwareObjects.mockup.BeamlineActionsMoc
    {"type": "annotated", "command": "HardwareObjects.mockup.BeamlineActionsMoc
  ]
</commands>
</object>
```

## Queue entry / task - For collecting data



Wavelength: 1.00 Å Detector: 10.1 mm Flux: 1.50e+10 ph/s UNKNOWN

Path: /data/id29/inhouse/opid291/20220916/RAW\_DATA/Mb/Mb-Mb/

Filename: Mb-Mb\_RUN#[IMG#]

Subdirectory: Mb/Mb-Mb/

Prefix: Mb-Mb

Acquisition

Exp Time: 0.02 Align Chip: ☒

Sub Sampling: 4 Take Pedestal: ☒

Default Parameters Run Now Add to Queue

Write a task that takes a Pydantic model and add it to available\_methods of Beamline object

```
82 legacy_parameters: LegacyParameters
83
84 class SsxChipCollectionQueueEntry(BaseQueueEntry):
85     """
86     Defines the behaviour of a data collection.
87     """
88     DATA_MODEL = SsxChipCollectionTaskParameters
89     NAME = "SSXChipCollection"
90     REQUIRES = ["point", "line", "no_shape", "chip", "mesh"]
91
92     # New style queue entry does not take view argument,
93     # adding kwargs for compatability, but they are unused
94     def __init__(self, data: SsxChipCollectionTaskParameters, view=None):
95         super().__init__(view=view, data_model=TaskNode(data))
96
97     def execute(self):
98         super().execute()
99
100
```

```
available_methods:
  datacollection: True
  characterisation: True
  helical: True
  xrf_spectrum: True
  energy_scan: True
  mesh: True
  ssx_chip_collection: True
  gphlworkflow: True
  test_collection: True
```




# NEW SAMPLES VIEW

MxCuBE-Web (OSC)

Get samples from SC    Filter:

Cell 1 Selected




Sample Changer


Cell 1 ☐ ☐ ☐

Puck 1	Puck 2	Puck 3
<input type="checkbox"/> Sample 1:1:01 1:1:01	<input type="checkbox"/> Sample 1:2:01 1:2:01	<input type="checkbox"/> Sample 1:3:01 1:3:01
<input type="checkbox"/> Sample 1:1:02 1:1:02	<input type="checkbox"/> Sample 1:2:02 1:2:02	<input type="checkbox"/> Sample 1:3:02 1:3:02
<input type="checkbox"/> Sample 1:1:03 1:1:03	<input type="checkbox"/> Sample 1:2:03 1:2:03	<input type="checkbox"/> Sample 1:3:03 1:3:03
<input type="checkbox"/> Sample 1:1:04 1:1:04	<input type="checkbox"/> Sample 1:2:04 1:2:04	<input type="checkbox"/> Sample 1:3:04 1:3:04
<input type="checkbox"/> Sample 1:1:05 1:1:05	<input type="checkbox"/> Sample 1:2:05 1:2:05	<input type="checkbox"/> Sample 1:3:05 1:3:05
<input type="checkbox"/> Sample 1:1:06 1:1:06	<input type="checkbox"/> Sample 1:2:06 1:2:06	<input type="checkbox"/> Sample 1:3:06 1:3:06
<input type="checkbox"/> Sample 1:1:07 1:1:07	<input type="checkbox"/> Sample 1:2:07 1:2:07	<input type="checkbox"/> Sample 1:3:07 1:3:07
<input type="checkbox"/> Sample 1:1:08 1:1:08	<input type="checkbox"/> Sample 1:2:08 1:2:08	<input type="checkbox"/> Sample 1:3:08 1:3:08
<input type="checkbox"/> Sample 1:1:09 1:1:09	<input type="checkbox"/> Sample 1:2:09 1:2:09	<input type="checkbox"/> Sample 1:3:09 1:3:09
<input type="checkbox"/> Sample 1:1:10 1:1:10	<input type="checkbox"/> Sample 1:2:10 1:2:10	<input type="checkbox"/> Sample 1:3:10 1:3:10
<input type="checkbox"/> Sample 1:1:11 1:1:11	<input type="checkbox"/> Sample 1:2:11 1:2:11	<input type="checkbox"/> Sample 1:3:11 1:3:11
<input type="checkbox"/> Sample 1:1:12 1:1:12	<input type="checkbox"/> Sample 1:2:12 1:2:12	<input type="checkbox"/> Sample 1:3:12 1:3:12



M Oscarsson and JB Florial




# IN SITU DATA COLLECTION




- Experimental workflow within the current implementation :
- Crystal screening:
  - based on known position X-Y CRIMS (synch) (Harvesting plan)
  - Or mesh scan on each drop







# BRAGGY INTERACTIVE MESH



Socket can be sent also to adxv,  
albulu


M Oscarsson and A Bocciarelli





**MXCuBE-W** **Chip**


Beamline Action

Beam size: 50






Chip (Diamond)  Navigate


Omega: 270.00  90.0

X: 0.177  0.1

Y: 15.657  0.1

Sample alignment

 Show motors

Chip

Move to

Add to queue

Equipment System log

Help Remote Sign out (opid291)

10.0 %

Detector UNKNOWN Sample Changer READY Capillary Exception [Errno 111] Connection refused Fast Shutter Exception [Errno 111] Connection refused Safety shutter CLOSED

Backlight Frontlight Video size

Run Queue Unmount Settings




Log messages:

- ✓ [10:59:06]: Not mounting next sample automatically (Auto mount next)
- ✓ [10:59:06]: Acquired (18.0, 15.8, -0.17, 22.0, 15.8, -0.17, 18.0, 24.0, -0.17)
- ✓ [10:59:06]: shutter closed
- ✓ [10:54:03]: Acquiring 49200



50  $\mu$ m

ESRF






# HIT-FINDING IN DAQ SEQUENCE



- **Identify** the presence of **diffraction spots** directly in the **DAQ**
- **LlmA2** detects spots and flags frames above a certain threshold directly in GPU
  - Possibility for automatic sorting reducing saved data and data bandwidth
  - Algorithm based on pyFAI is the same as in the offline analysis
  - Currently under final commissioning
  - Hits are written and flagged in the hdf5 files directly


S. Debionne, A. Homs Puron, J. Keiffer, R. Ponsard, L. Claustre, A. Götz, P. Fajardo  
European Synchrotron Radiation Facility (ESRF), Grenoble, France

- LImA2: evolution of LImA towards distributed systems (from scratch)
- Goals: **scalable** DAQ and **low latency** processing



- C++ library
- Separate DAQ from processing
- **Feature extraction** and **advanced data reduction**
- Specific & optimized pipelines






Parallel streams:

- 2x computers
- 2x receivers / computer


GPU:

- Geometry assembly
- Pedestal/gain correction
- Data sparsification
- Peak finding

- Sparse data: x4-x100 factor reduction
- *Is-hit* information can be used as veto
- In production since January 2023
- Low-latency online visualization
- Continuous processing @ 500 Hz  $\Rightarrow$  1 kHz



- Migrate to int16 or float16
- Reject frames on the fly
  - Reduce file sizes and increase frame rate
- Off-line image sorting
  - Keep only meaningful frames in the long term



- Advance on SSX ID29 specific features
- Develop specific data sorting / pre-analysis
- Braggy visualisation of Jungfrau images
  - Hits display
  - Max Projection
- Deploy MASSIF1 workflows on other beamlines
- Port data collection routines to Bliss
- ICAT integration
- Deploy SSO (and 2FA) and proposal selection